



BRAVO

Equipos médicos y terapéuticos





10D Rotating Green
Laser Lights

Newest Painless Fat
Removal Choice

From **MaxMaster Slim**



Fat loss in a healthier way

»» As the indispensable part of human's endocrine system, fat cells play a critical role in sensing and responding to changes in systemic energy balance and the metabolic cycle. MaxMaster Slim does not kill them like liposuction and cryolipolysis. It may sound reasonable to kill the fat cells in weight loss but actually it is harmful and has side-effects. When the fat cells are destroyed, the future fat will be stored and accumulated in other area and may cause asymmetric accumulation of fat.

»» MaxMaster Slim opens a new possibility in cellulite reduction.

By using the low level laser, fat cells can be shrinked and then consumed as energy and flushed out naturally via the body's lymphatic system when combined with light exercise and sufficient water intake. It is non-invasive and no side-effects reported.




Low Level Laser Therapy for Fat Loss

Safest




Comfortable



No Downtime



No Freezing




No Bruising



No Swelling

Fast results (most patients see results in 4 weeks or less)



10 Dancing Laser Light

360° Scanning on Larger Area

Shorten Treatment Time and
Improve the Fat Reduction Efficiency



Upgraded Laser Head

-Multi diode laser foldable on both wings can focus multiple lasers and make the slim treatment more flexible

High Quality Diode Laser

- Imported from Japan
- Double lens to increase the spot area of the laser and concentrate the energy
- High-temperature-resistant materials for continuous running

Industrial-grade Elevator

-Make the entire laser head fully automatically up and down to adjust to the height of the bed

12.1 Inch high-sensitivity Touch Screen

- User-friendly
- More Convenient




About the Treatment



MaxMaster Slim provides a quick and no downtime treatment. We can resume normal daily activities immediately after the treatment.

We just need to simply lie under the laser on the front for (20) minutes and then on the back for (20) minutes.

Easy!!





TECHNICAL PARAMETERS

Maximum reach of laser head	110cm
Angle adjustable of laser wings	100 degree
Weight of laser head	12kg
Maximum reach of Elevator	500mm
Size of screen	12.1 inches
Power of diode	200mw
Wavelength of diode	532nm
Voltage	90v-240v
Power	120w
Number of diode	10pcs

10D

LOW LEVEL LASER THERAPY

For Fat Loss

